


# Total Lift Bed



## TILTING INSTRUCTIONS AND TRACKER

1. Ensure clearance of all lines/leads. Appropriate socks/footwear as indicated.
2. The strapping system is to be firmly fastened during tilting to ensure patient safety. Strap placement is based on patient's clinical status. For example, straps can be placed just above and below the knees, at the chest or directly over the knees and at the pelvis and chest. Ask a super user for alternate strap configurations if needed.
3. For tilting, hold down the bottom two buttons (**11 and 12**) simultaneously until beep is heard to unlock tilt mode.
4. Then hold down bottom right button to begin tilting (**12**).
5. As the footplate moves in to meet the patient, ensure proper positioning of the patient's feet/knees.
6. The bed will then go flat (on newer models it does not have to go completely flat) and slowly increase in tilt position. Be sure to watch for clearance of all lines for safety. The inclinometer is on the right side of the bed under the top bed rail and/or on the upper siderails.
7. Tilt to the angle patient last tolerated, allow patient to acclimate and then increase angle if able. Mark the time the tilt session began and write the total amount of time on the chart below.
8. If the patient is tolerating well, increase the indicated degree based on the below listed parameters.
9. When the session is complete, hold down the CPR flat button until the patient returns to zero degrees.
10. You can now position the patient as warranted. Ensure mattress low air-loss and alternating pressure features are still correct and patient appropriate.



**MANUAL  
EMERGENCY  
LOWERING**

LOCATED PATIENT RIGHT SIDE,  
HEAD OF BED, BEHIND SCALE  
BOX. REMOVE PIN AND PULL  
LEVER.

DATE	START TIME	STOP TIME	DEGREE	WEIGHT BEARING	INITIAL



TILT INSTRUCTIONS VIDEO



HAND PENDANT INSTRUCTIONS



MORE TRAINING RESOURCES

